

فناوران نانو مقیاس



ماسک های تنفسی نانو


نادر نادری

مدیر عامل شرکت فناوران نانو مقیاس



پنجشنبه ۲۹ آبان ماه ۱۳۹۹

Nanofibers



4.34 nanofibre

nano-object with two similar external dimensions in the nanoscale and the third dimension significantly larger


NOTE 1 A nanofibre can be flexible or rigid.

NOTE 2 The two similar external dimensions are considered to differ in size by less than three times and the significantly larger external dimension is considered to differ from the other two by more than three times.


NOTE 3 The largest external dimension is not necessarily in the nanoscale.


[ISO TS 27687:2008, definition 4.3]

Electrospinning



The velocity of the jet can reach values of a few m/s, and strain rates are up to 10^7 s^{-1} .





About us



Fanavaran Nano-Meghyas (FNM Co. Ltd.;) was founded in 2004, is a knowledge based company and its goals are the development of nanofibers technology and its applications. FNM's products and services are design and production of electrospinning machines in lab, pilot and industrial scales as well as melt, force and blown electrospinning systems, with various accessories (High Voltage power supplies, Syringe Pumps and collectors), with focus on producing of respiratory face mask, power plant and automotive air and oil filters, window screen, vacuum cleaner bags, cosmetic face masks, wound dressing and etc. based on electrospun nanofibers.

Name of Company:	Fanavaran Nano-Meghyas (Fnm co. Ltd.)
E-mail:	fnm.ir.co@gmail.com, info@fnm.ir
Office:	Unit 10, No. 23, 2nd Keyhan Ave., Ayatollah Kashani St., Tehran, Iran,
Factory:	West 4th Street, Golgoon Industrial zone, 5th km Shahriar Road, Karaj old Road. Iran
Nanofiber Innovation Center	No. 8; Hamedan Ave. North Kargar Street; Tehran; Iran
Tel.:	+98 (21) 65612497
ICANANO	ICAN; NANO site; Sh. Ehsani Rad ST., Engelab St., Parsa Sq., Ahmadabad Mostoufi Rd., Azadegan Highway, Tehran, Iran


Awards, Verification and Standards

- Ranked eighth in International Nanotechnology Festival-Tehran (2009).
- Award for research project at the Festival of Science to Practice (December 2010)
- Third laureate R&D 25th Khwarazmi international Award, 5th Feb 2012 Tehran, Iran.
- Academy of science award in developing countries (TWAS); UNESCO, 2012
- Second rank in Technology at 6th National Nano-Awards Festival, October, 2011, Tehran, Iran.
- First laureate nano products award at 11th National Nano-Awards Festival, October, 2016, Tehran, Iran.
- Award for high tech export at 12th National Nano-Awards Festival, October, 2017, Tehran, Iran.
- Award for high tech export at 13th National Nano-Awards Festival, November, 2018, Tehran, Iran.
- 5 International and 12 Iranian Patents





Today there are **8 overseas offices/sale representative**, in Shanghai (China), Suzhou (China), Kuala-Lumpur (Malaysia), Jakarta (Indonesia), London (United Kingdom), Istanbul (Turkey), and Islamabad (Pakistan).




Citations to our Products in Research Papers




More than 250 published research papers (updated 2018).

Citations are available on our website

(<http://en.fnm.ir/papers/products-papers>)



Number of publications using FNM Co. products (updated: August 2020)



Company Facilities



Analysis Equipment

- Scanning Electron Microscopy (SEM)



- FNM Filter Test machine



FNM Filter/Mask Test Machine



Specifications

		FT150EA	FT200PO	FT200PS
Test Modes	Pressure drop test	✓	✓	✓
	Air permeability test	✓	✓	✓
	Filtration Efficiency Test	Atmospheric	Oil Particles	Oil and Salt Particles
	Bubble Point	×	Optional	Optional
Standards	BS EN 149	✓	✓	✓
	BS EN 779	✓	✓	✓
	ISO 16890	✓	✓	✓
	ISO 16900-3	✓	✓	✓
	ISO 11155-1	✓	✓	✓
	ISO 5011	✓	✓	✓
Control	PLC	✓	✓	✓
	HMI	7"	7"	7"
	Ambient Temp.	×	✓	✓
Air flow	Flow	10 - 150 l/min	10 - 200 l/min	10 - 200 l/min
	Digital control	✓	✓	✓
Media Holder	Area	100 cm ²	20, 50 and 100 cm ²	20, 50 and 100 cm ²
Sensors	Temperature	✓	✓	✓
	Relative Humidity	✓	✓	✓
	Digital Tank Pressure Control	×	✓	✓
	Tank Pressure	1 - 8 bar	1 - 8 bar	1 - 8 bar
Pressure Drop	Pressure Drop	0 - 1200 Pa	0 - 1200 Pa	0 - 1200 Pa
	Digital control	✓	✓	✓
Air Dryer	Air Line Trap	✓	✓	✓
	Dryer System	×	✓	✓
Particle Counter	Laser Particle Counter	1	2	2
	Channels	six-channel	six-channel	six-channel
	Channel Sizes	0.3, 0.5, 1, 2.5, 5, 10 μm	0.3, 0.5, 1, 2.5, 5, 10 μm	0.3, 0.5, 1, 2.5, 5, 10 μm
	Flow rate	2.8 L/min	2.8 L/min	2.8 L/min
Aerosol Generator (Oil)	Generator	×	✓	✓
Aerosol Generator (Salt)	Generator	×	×	✓
	Neutralizer	×	×	✓
Aerosol Dilutor	Particle Dilution	×	100:1	100:1
	Type of Aerosol Challenge	×	PSL, PAO, DOP	PSL, PAO, DOP, NaCl
Printer		Optional	Optional	✓
Respiratory Face Mask Holder		Optional	Optional	Optional
Power		Single phase, 220 V, AC	Single phase, 220 V, AC	Single phase, 220 V, AC
Weight (kg)		About 170 kg		
Size, (Length, Width, Height)		94 cm, 92 cm, 163 cm		

BS EN 149:2001 +A1:2009



DATE: 13 February, 2020
NO.: 98112482

Respiratory Face Mask Filtration performance test results


Requested by: Zist Abzar Pajooohan Co.

Test Conditions

- Test apparatus: FNM Filter Test
- Test aerosol: Oil
- Air flow rate for efficiency test: 32 l/min
- Sample name: **FFP2 - Rima**
- Particles tested size: 0.3, 0.5, 1.0, 2.5, 5.0 and 10.0 μm

Filtration performance of respiratory face mask

Performance results (According to BS EN 149)								
NO.	Efficiency (%)						Pressure Drop (Pa) @ (l/min)	
	0.3 μm	0.5 μm	1.0 μm	2.5 μm	5.0 μm	10.0 μm	30	95
1	96	97	98	99	100	100	52	187
2	95	95	96	99	100	100	51	185
3	96	98	99	100	100	100	52	187



DATE: 13 February, 2020
NO.: 9811241

Respiratory Face Mask Filtration performance test results


Requested by: Zist Abzar Pajooohan Co.

Test Conditions


- Test apparatus: FNM Filter Test
- Test aerosol: Oil
- Air flow rate for efficiency test: 32 l/min
- Sample name: **FFP3 - Rima**
- Particles tested size: 0.3, 0.5, 1.0, 2.5, 5.0 and 10.0 μm

Filtration performance of respiratory face mask

Performance results (According to BS EN 149)								
NO.	Efficiency (%)						Pressure Drop (Pa) @ (l/min)	
	0.3 μm	0.5 μm	1.0 μm	2.5 μm	5.0 μm	10.0 μm	30	95
1	99	99	100	100	100	100	63	231
2	99	99	99	100	100	100	64	228
3	99	99	99	99	100	100	63	230



<http://icanano.ir/>




پلتفرم نانوالیاف

New Standard in ISO

Nanotechnologies –Air filter media containing polymeric nanofibres;
Specification of characteristics and measurement methods

- 1- Potential New Work Item Proposal: **WG4 Interim meeting in Korea, June 2015**
- 2- Circulation of NWIP: **2015, June, 30**
- 3- Result of Ballot (approved): **2015, Sept., 24**
- 4- Interim Meeting: **May 2016 Kyoto**
- 5- Interim Meeting: **November 2016 Singapore**
- 6- web meeting: **May 2017**
- 7- Interim Meeting: **November 2017 Seoul**
- 8- web meeting: **Feb. 2018 (resolve the USA comments)**
- 9- web meeting: **July 2018 (resolve the USA comments)**
- 10- Interim Meeting: **November 2018 Malaysia**
- 11- new version was circulated for DTS: **Feb. 2019**
- 12- Result of DTS Ballot: **approved 2019, April, 16**
- 13- Meeting Sydney, **May, 2019 for resolving comments**
- 14- 4 web meeting: **July 2019 (resolve the Japan comments)**

ISO/TS 21237




Lab-Scale Electrospinning unit



Hybrid Electrospinning machine

۱۴



**Cartridge and needle-based
Electrospinning machine**

Pilot-scale Electrospinning Machine

۱۵



Industrial Electrosinning Machine

۱۶




Yield



Up to 100 g/hr


**depending on polymer solution and electrospinning parameters*

Nanofiber sizes of
60-500 nm*



BATH/DIP ELECTROSPINNING






CONVENTIONAL MULTINEEDLE ELECTROSPINNING




http://www.bjtechnova.com/product/show.asp?lb_id=8
<https://www.qingzitech.net/electrospinning-nanofiber-production-line-mf01>

Upscaling: what to think about?



- **Production rate**
- **Wastage**
- **Maintenance/Cleaning**
- **Flexibility**
- **Adhesion to support layer**
- **Solution Concentration**
- **Flammability of solvents**
- **Thick nanofiber mat**

Blowing-Assisted Electrospinning



Comparison of different electrospinning methods

Method	Production rate*	CONSTANT Solution Concentration	Flexibility (variety of materials)	Solution waste**	flammability	Adhesion to support layer	Thick nanofiber mat
Multineedle electrospinning	Poor 2gr/h/unit (10 needles)	Excellent	Good	Very Low <5%	Very Low	Low	Medium
Bath/Dip electrospinning	Medium 3gr/h/unit	Poor	Poor	High >50%	High	low	Low
Cartridge electrospinning	Medium 4gr/h/unit	Poor	Medium	Medium >25%	Medium	Low	Low
Blowing- assisted electrospinning	High 12gr/h/unit	Excellent	Excellent	Very low <5%	Very Low	High	High Up to 2 mm

*Polyamide

**Based on repeated internal testing

Filter Application


- HVAC air filter
- Gas turbine air filter
- Automotive air filter
- Clean room air filter
- Bag house air filter
- Antimicrobial air filter
- Catalytic filter
- Highly selective filter
- Respiratory face masks




Filter Application

Major Performance Characteristics for Air Filters based on Market Survey

- 1-Efficiency
- 2-Pressure Drop
- 3- Dust holding Capacity (Life Time)





Microfiber of main filters paper



Filters' paper classification



Class	Final Pressure Drop Pa	Average efficiency (Em) of 0,4 μm particles %
F7	450	$80 \leq Em < 90$
F8	450	$90 \leq Em < 95$
F9	450	$95 \leq Em$

Filtration




× 40

SEM Images of Nanofibers on filter paper



EMCRAFTS	WD	HV	Probe	1.0 μm	x10000	2016.10.26	
	5.37mm	20.00kV	1.00			DET	SEI





FNM	WD	HV	Probe	10.0 μm	x975	2018.7.4	
	26.72mm	20.00kV	4.00			DET	SEI



Filtration Efficiency tests

Active nanofibrous filtering layer with narrow pore size suitable for submicron filtration

A layer of NF changes the



EN 779:2002. AIR FILTER TEST RESULTS

GENERAL

Test no.:	122876	Date of test:	28.6 and 16.- 17.7.2012	Supervisor:	RHo
Test requested by:	Behran Filter Co.			Device receiving date	16.5.2012
Device delivered by:	Behran Filter Co.				

DEVICE TESTED

Model	Manufacturer	Construction
Gas Turbine Air Filter V94.2	Behran Filter Co.	Cylindrical filter
Type of media	Net effective filtering area	Filter dimensions (diameter x length) (The length includes gaskets)
90/40 EPE K WB2-G+NANO	19 m ²	328 mm x 624 mm

TEST DATA

Test air flow rate	Test air temperature	Test air relative humidity	Test aerosol	Loading dust
0.347 m ³ /s	24 - 25 °C	38 - 43 %	DEHS	ASHRAE

RESULTS


Initial pressure drop	Initial arrestance	Initial efficiency (0,4 µm)	Dust holding capacity	Untreated / discharged efficiency of filter material (0,4 µm)
238 Pa	>99 %	72 %	8 / 226 / 504 g	Non Applicable
Final pressure drop	Average arrestance	Average efficiency (0,4 µm)	Filter class (450 Pa)	
250 / 350 / 450 Pa	>99 / >99 / >99 %	79±1 / 97±0 / 99±0 %	F9 (0.347 m ³ /s)	

Remarks: -

NOTE: The performance results cannot by themselves be quantitatively applied to predict filter performance in service.
The results relate only to the tested item.

THE **INVISIBLE KILLER**

Air pollution may not always be visible, but it can be deadly.





Air pollution was estimated to cause 4.2 million premature deaths worldwide per year.
(April of 2019, www.who.int/airpollution/en/)






Respiratory Face Mask



Respiratory Face Mask



Nanofiber-based Anti-Allergy Bedding Test Results




SEM image of nanofiber-based anti-allergy bedding




Nanofiber-based anti-allergy bedding substrate

Performance of Nanofiber-based Anti-Allergy Bedding

Performance of nanofiber-based anti-allergy bedding (According to EN 779, EN ISO 9237)							
Sample Name	Efficiency (%)						Air Permeability (l/m ² /s)
	0.3 μm	0.5 μm	1.0 μm	2.5 μm	5.0 μm	10.0 μm	
Nanofiber-based Anti-Allergy Bedding	72	77	81	88	100	100	181 (@ 125 Pa)
PP Spunbond Substrate	11	12	35	45	50	56	250 (@ 60 Pa)



Nanofiber-based Beauty Face Mask Test Results



SEM image of nanofiber-based
beauty face mask




Beauty Face Mask

Release Test of Nanofiber-based Beauty Face Mask






Test Method	λ_{\max}
Spectrophotometry	255 nm



Lithium Ion Battery Separator

FNM Separators in Lithium-Ion Battery

1. The porosity is more than 60%
2. Electrolyte uptake 500%
3. Thermal stability above 180 ° C
4. High ionic conductance
5. No battery electrolyte leakage
6. Higher safety offered by high heat resistance
7. Shorter processing time for electrolyte injection




Sample	Electrolyte uptake (%)	Porosity (%)	Ionic conductivity (mS.cm ⁻¹)
PP	86	36	0.8
FNM	500	65.02	2.3




Lithium Ion Battery Separator

Impedance




Electrochemical stability



Cycle Performance

The cycle performance of the cells using FNM separator is studied at discharge rate of 0.5C as a function of cycle number (up to 200 cycles).



Behran Filter;
4 Industrial and 1 pilot electrospinning machines

Azad Filter; (Industrial filters)
1 Industrial electrospinning machine

Serkan Filter (Car filters)
1 Industrial and 1 pilot electrospinning machines

Zist abzar Co. (Respiratory face masks)
2 Industrial and 1 pilot electrospinning machines

Deylaman Filter (Industrial filters)
1 Industrial electrospinning machines

Golrang Group (Respiratory face masks)
1 Industrial electrospinning machines






**Malaysia: 3 Pilot
electrospinning
machines, 20 lab-scale**

*A Joint Company in Malaysia:
NanoLab Instruments*




South Korea; 1 Industrial electrospinning machine



China: 6 Industrial, 5 lab-scale electrospinning machines and 1 pilot electrospinning machines



صادرات دستگاه تست فیلتر به کشور تایلند



اطلاعات تماس شرکت

موبایل: ۰۹۱۲۶۲۷۱۸۲۸

ایمیل: fnm.ir.co@gmail.com

تلفن: ۶۵۶۱۲۴۹۶


وبسایت شرکت: fnm.ir

اینستاگرام: [rimamask](https://www.instagram.com/rimamask)

تلگرام:

واتساپ: ۰۹۱۲۶۲۷۱۸۲۸

تشکر از توجه شما



FNM Co.

en.fnm.ir

